

Technical Sheet

The ALN Extraction and/or Repositioning Kits


- Only by jugular approach
- Pincer system: specially designed for the ALN Optional Vena Cava Filter
- CE Mark 0297
- FDA Cleared

ALN Straight RS & ALN Pre-curved RS


CLASS III (according to Directive 93/42/EEC)

The ALN Extraction and/or repositioning kits ALN STRAIGHT RS and ALN PRE-CURVED RS contain:

- ✓ A **Pincer catheter** (7F sheath with two black markers + pincer with 8 strands + 3-ways stop cock)
- ✓ An **introduction system** (dilator, 9F extraction sheath with a radio-opaque marker band)
- ✓ A **puncture needle**
- ✓ A **J-guidewire**










Parameters of the Retrieval System

Diameter of the strands (mm)	Material	Diameter of the pincer (mm)	Length of the pincer (mm)
0.4	Stainless 316 LVM*	13.5-15	24 








*ISO 5832-1

The ALN Extraction and/or Repositioning Kit ALN Straight RS

ALN STRAIGHT RS	
Approach	Jugular
Product Code	FT.902010
Introduction system: Dilator in the extraction sheath	
	
Extraction sheath PEHD Compound PEHD	 Radio-opaque marker band (Bismuth) Ø Internal: 9F Length: 570 mm No haemostatic valve Ø External: 11F
Dilator PEHD Compound ABS / PEHD	 Ø Internal: 3F Length: 630 mm Ø External: 9F
Sheath PEHD Compound PEHD	 Two black markers to visualize the stages of the removal procedure Ø Internal: 7F Length: 596 mm Ø External: 9F
Pincer Kit	 Length: 730 mm 7F Sheath internal diameter Handle 3-ways stop cock Pincer
Puncture needle Stainless 304 PE / TPX	 Ø : 17 G
J-guidewire Stainless 304 with teflon coating	Ø : 0.035" Ø Tip: 9F (flexible J) Length: 150 cm 

Sterilization by ethylene oxide

The ALN Extraction and/or Repositioning Kit ALN Pre-curved RS

ALN PRE-CURVED RS	
Approach	Jugular
Product Code	FT.902010/VS2
<p>Introduction system: Dilator in the extraction sheath</p> 	
Extraction sheath PEHD Compound PEHD	 <p>Radio-opaque marker band (Bismuth) No haemostatic valve</p> <p>Ø Internal: 9F Ø External: 11F</p> <p>Length: 570 mm</p>
Dilator PEHD Compound ABS / PEHD	 <p>Ø Internal: 3F Ø External: 9F</p> <p>Length: 630 mm</p>
Sheath PEHD Compound PEHD	 <p>Two black markers to visualize the stages of the removal procedure</p> <p>Ø Internal : 7F Ø External : 9F</p> <p>Length: 596 mm</p>
Pincer Kit	 <p>Length: 730 mm</p> <p>Handle 3-ways stop cock 7F Sheath internal diameter Pincer</p>
Puncture needle Stainless 304 PE / TPX	 <p>Ø : 17 G</p>
J-guidewire Stainless 304 with teflon coating	<p>Ø : 0.035" Ø Tip: 9F (flexible J) Length: 150 cm</p> 

Sterilization by ethylene oxide